


INTRODUCTION



Phage production's downstream process includes a critical purification step susceptible to destabilise phage suspensions. Ideally, phage titer should be determined before fill and finish to avoid production losses, but titration is a time-consuming method. The colloidal properties of phages can be used to rapidly determine the concentration and size distribution of suspended particles, e.g., with interferometric light microscopy (ILM) method. **Can ILM successfully monitor phage purification ?**

MATERIALS AND METHODS

Phage purification

1 Frontal filtration

Purification intermediate (PI)


2 Tangential flow filtration

Diluted PI (dPI)

Washed PI (wPI)

Concentrated PI (cPI)

Formulated PI (fPI)




Phage anti
Staphylococcus aureus
vB_SauM-V1SA20
Family: *Herelleviridae*
Genus: *Silviavirus*
Myovirus

N=3 purifications

50 μ L sampled at each steps:

- Phage titer (N=3 spot –test measures)
- Particle concentration/titer and size (N=2 ILM measures)

Interferometric light microscopy (ILM)



Camera detects interference patterns



No treatment





7 μ L sample

White light (LED)

Application range:

Concentration: $10^8 - 10^{10}$ part/mL
Particle size: 80 nm – 1 μ m

RESULTS



● PI ○ dPI □ wPI △ cPI × fPI — Linear regression

Deviation = $(\log(\text{particle concentration}) - \log(\text{titer})) / \log(\text{titer})$

Deviation	PI	dPI	wPI	cPI	fPI
Purification 1	5 %	4 %	5 %	3 %	3 %
Purification 2	4 %	6 %	7 %	4 %	4 %
Purification 3	4 %	6 %	6 %	10 %	4 %
Mean $\pm$ SD	$(4 \pm 1) \%$	$(6 \pm 1) \%$	$(6 \pm 1) \%$	$(6 \pm 4) \%$	$(4 \pm <1) \%$

CONCLUSION

Assets:

Small sample (5-10 μ L) At any viscosity
No necessary treatment Rapid measure (a few minutes)

Monitoring of purification:

- ✓ Does the suspension become homogeneous?
- ✓ Is the physical titer preserved?
- ✓ Does the size distribution narrow?

Linear correlation $\Rightarrow$ physical titer may estimate phage titer

Limits:

Highly concentrated samples ($>10^8$ part/mL)
The higher the concentration, the more precise the result (lowest deviation of 4 % reached for PI and fPI)
Detects particle > 80 nm $\Rightarrow$ **not for podovirus** morphotype
Measures particles, not phages:
 $\Rightarrow$ **Cannot detect phages with emptied capsid**
 $\Rightarrow$ **Cannot replace titration as release quality control**

ACKNOWLEDGEMENTS